



Zenith 11 Medical Catalog

Elevating Healthcare Excellence



Company Overview

Zenith 11 Medical: Elevating Healthcare Excellence

At Zenith 11 Medical, we are committed to enhancing the quality and standards of healthcare services. Our dedication to reaching new heights in medical care is reflected in our embrace of innovation, advanced technologies, and continuous improvement. We strive to ensure superior patient outcomes through a holistic approach to healthcare that prioritizes excellence in medical practice, patient experience, and overall wellness. Whether through cutting-edge treatments, compassionate care, or the pursuit of medical advancements, we uphold the highest standards in the healthcare industry.





Our Mission

At Zenith 11 Medical Device Company, our mission is to revolutionize healthcare through innovation. We are committed to designing and delivering state-of-the-art medical devices that enhance patient care and outcomes. Our focus is on developing accessible, efficient, and reliable solutions that meet the evolving needs of healthcare professionals and patients alike. By combining advanced technology with user-friendly designs, we strive to make quality healthcare more attainable and effective for communities worldwide.

Zenith 11 Medical

Products & Services

PUSH CLICK Y CONNECTOR KIT

Large inner cavity passage diameter to meet the needs of various instrument delivery

Category: Push Click Y Connector Kit




FEATURES

- 01 Special switch design can automatically rebound to close, effectively reducing blood loss
- 02 Enlarged knob, threaded design, more comfortable pressing
- 03 Easy one-handed operation Maximum pressure is 180psi
- 04 The maximum inner diameter is $2.7\pm0.1\text{mm}$, which can meet the delivery of more instruments Introducer and Torque for 0.014"-0.021" guide wires

Description	Ref. No.	Torque Device	Guidewire Introducer	Extension Tubing
Push-Click Y Connector	DPY01	N	N	N
	DPY01	N	N	Y
Push-Click Y Connector	DPY01	Y	Y	N
	DPY01	Y	Y	Y

Products & Services



RADIAL ARTERY COMPRESSION DEVICE

Easy Wearing and Stable Pressure
Control

Category: Radial Artery Compression Device



FEATURES

- 01 Transparent plate for easy observation and precise location on the puncture point
- 02 Soft elastic material for patient's comfort
- 03 Force concentration by spherical design for effective hemostasis
- 04 Non-touch with wrist to keep from harm on other vessels
- 05 Self-lock Fixed Loop for Easy Wearing and Stable Control on Pressure
- 06 Helpful Elastic Band with Scale for Releasing Gradually

Description	Ref. No.	Length of Plate (mm)	Length of Band (mm)
	DPY01	90	140
	DPY01	90	140

Products & Services

BALLOON IN-DEFLATION DEVICE

Patent Selok Self-Locking Mechanism with Maximum Pressure of 30atm

Category: Balloon In-deflation Device



FEATURES

01 Surface appearance patent and ergonomic design

02 Selok Patent technology

03 20ml clear syringe barrel

04 Precise pressure gauges



BENEFITS

- Suitable for the left and right hand, easy to operate, reduce fatigue
- Stable pressure output
- Fast and easy pressure release
- High pressure display accuracy

Description	Ref. No.	Nominal Max Pressure (atm)	Capacity(ml)	Include 3-Way Stopcock
Mastro Balloon in-Deflation Device	DID30S(OFF)	30	20	YES

Products & Services

ARM SUPPORT BOARD

Carbon fibre, unibody, lightweight and durable
Good development effect under X-ray
Ergonomically designed to conform to the patient's arm

Category: Arm Support Board



FEATURES

- 01 Carbon fibre, unibody, lightweight and durable
- 02 Good development effect under X-ray
- 03 Ergonomically designed to conform to the patient's arm


Ref.No.	Size (mm)
DBT03R	670*580
DFJ01	365*595
DBT03L	670*580

Products & Services

3-PORT MANIFOLD

3-Port Manifold with Pressure Rating Up to 500psi Clear Body Allowing Easy Visibility for De-bubbling

Category: Manifolds



FEATURES

- 01 3-Port Manifold with Pressure Rating Up to 500psi
- 02 Clear Body Allowing Easy Visibility for De-bubbling
- 03 Easy Rotating Handles
- 04 Rotating Luer for Secure Connections
- 05 Optional Extension Tubing
- 06 ON/OFF Handle for Options

Description	Ref. No.	Length of Plate (mm)	Type of Stopcock
3-Port Manifold	DM5003(ON)	500psi	ON
	DM5003(OFF)	500psi	OFF

Products & Services

3-WAY STOPCOCK

Carbon fibre, unibody, lightweight and durable
Good development effect under X-ray
Ergonomically designed to conform to the patient's arm

Category: Manifolds



FEATURES

- 01 With a Male Luer Lock and Two Female Luer Connectors
- 02 Pressure Rating Up to 500psi
- 03 Poly-Carbonate Clear Body for Easy Visualization of Fluid Flow
- 04 Smooth Turning Handle with “ON” or “OFF” Options

Description	Ref. No.	Pressure Resistance	Type of Stopcock
3-Way of Stopcock cock	DM5003(ON)	500psi	ON
	DM5003(OFF)	500psi	OFF

Products & Services

Y CONNECTOR

Easy Rotating Handles

Category: Twist Y Connector Pack




FEATURES

- 01 Easy Rotating Handles
- 02 More options, Y-shaped, and T-shaped side branches, to meet different surgical needs

Name	Ref. No.	Torque device	Guidewire Introducer	Extension Tubing	Pressure(psi)	Shape
Y Connector Kit (Twist)	DHRY01B	Blue Torque	N	N	500	Y
Y Connector Kit (Short-Twist)	DHDY01B	Blue Torque	N	N	500	Y
Y Connector Kit (Short-Twist)	DHDY01B	Blue Torque	N	N	500	Y


Products & Services



INTRODUCER SHEATH

Enhanced Pushability and Flexibility
Enlarged Inner Diameter for Ultra-smooth
Transition Between Dilator and Sheath

Category: Introducer Sheath



FEATURES

- 01

Enhanced Pushability and Flexibility
- 01


Smooth transition section to effectively reduce damage to blood vessels
- 01

Specially Designed for Vascular Applications
- 01

Enlarged Inner Diameter for Ultra-smooth Transition Between Dilator and Sheath

Introducer Sheath Kit (Radial/Femoral)					
Description	Ref. No.	French Size	Sheath length (cm)	Guidewire Style	Needle Gauge
Radial Introducer Sheath	DQ05082121T/S	5F	8	0.021"	21
	DQ05112121T/S	5F	11		
	DQ06082121T/S	6F	8		
	DQ06112121T/S	6F	11		
	DQ07082121T/S	7F	8		
	DQ07112121T/S	7F	11		
Femoral Introducer Sheath	DQ05113818T	5F	11	0.038"	18
	DQ06113818T	6F			
	DQ07113818T	7F			
	DQ08113818T	8F			


Products & Services



INTRODUCER SHEATH (DXQ)

Polymer Guidewire with Hydrophilic Coating

Category: Introducer Sheath(Hydrophilic)



FEATURES

- 01

Enhanced Pushability and Flexibility
- 01

Smooth transition section to effectively reduce damage to blood vessels
- 01

Specially designed Valve to prevent blood flow back and air entry
- 01

Enlarged Inner Diameter for Ultra-smooth Transition Between Dilator and Sheath

Ref. No.	French Size	Sheath length (cm)	Guidewire Style	Needle Gauge
DXQ05082520	5F	8	0.025"	20G-30
DXQ05112520	5F	11		
DXQ06082520	6F	8		
DXQ06112520	6F	11		
DXQ07082520	7F	8		
DXQ07112520	7F	11		
DXQ08082520	8F	8		
DXQ08112520	8F	11		

Contact Information



510-828-8704



parth@zenith11.com



<https://zenith11medical.com/>



99198 Overseas Highway, Ste #12, Key Largo, Florida - 33037

